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1. Lloyd-Jones D, Adams RJ, Brown TM, et al. Heart Disease and Stroke Statistics—2010 Update. A Report from the American Heart Association Statistics Committee and Stroke Statistics Subcommittee. Circulation. 2010;121:e1-e170.

2. Wu A. B-Type natriuretic peptide and its clinical utility in patients with heart failure. *Medical Laboratory Observer* 2001; 10: 10-14.
3. Bonow RO. New insights into the cardiac natriuretic peptides. *Circulation* 1996; 93: 1946-1950.
4. McDowell G, Shaw C, Buchanan K, Nicholls D. The natriuretic peptide family. *Eur. J. Clin. Invest.* 1995; 25:291-298.
5. Yandle T. Biochemistry of natriuretic peptides. *J. Internal Med.* 1994; 235:561-576.
6. Mukoyama M, Nakao K, Hosoda K, Suga S, Saito Y, Ogawa Y, Shirakami G, Jougaski M, Obata K, Yasue H, Kambayashi Y, Inouye K, Imura H. Brain natriuretic peptide as a novel cardiac hormone in humans: Evidence for an exquisite dual natriuretic peptide system, atrial natriuretic peptide and brain natriuretic peptide. *J. Clin Invest.* 1991; 87:1402-1412.
7. Clerico A, Iervasi G, Del Chicca MG, Emdin M, Maffei S, Nannipieri M, Sabatino L, Forini F, Manfredi C, Donato L. Circulating levels of cardiac natriuretic peptides (ANP and BNP) measured by highly sensitive and specific immunoradiometric assays in normal subjects and in patients with different degrees of heart failure. *J. Endocrinol. Invest.* 1998; 21:170-179.
8. deLemos JA, Morrow DA, Bentley JH, Omland T, Sabatine MS, McCabe CH, H 全部 C, Cannon CP, Braunwald E. The prognostic value of B-type natriuretic peptide in patients with acute coronary syndromes. *New Engl. J. Med.* 2001; 345: 1014-1021.
9. Manual Guide – Safety Man 年龄 ment, No. CDC-22, Decontamination of Laboratory Sink Drains to Remove Azide Salts. April 30, 1976. Atlanta GA: Centers for Disease Control.
10. Approved Guideline – Procedures for Handling and Processing of Blood Specimens, H18-A2. 1999. National Committee for Clinical Standards.
11. Cembrowski GS, Carey RN. Laboratory Quality Man 年龄 ment: QC & QA. ASCP Press, Chicago, IL, 1989.
12. Hansen HJ, et al. HAMA interference with murine monoclonal antibody-based immunoassays. *J Clin Immunoassay* 1993; 16: 294-299.
13. Levinson SS. The nature of heterophilic antibodies and the role in immunoassay interference. *J Clin Immunoassay* 1992; 15: 108-115.
14. Approved Guideline Evaluation of precision performance of clinical chemistry device, EP5-A. 1999. National Committee for Clinical Laboratory Standards, V19, N2.
15. Krouwer JS, Rabinowitz R. How to improve estimates of imprecision. *Clinical Chemistry* 1984; 30: 290-292.
16. Maeda K, Tsutomoto T, Wada A, Hisanaga T, Kinoshita M. Plasma brain natriuretic peptide as a biochemical marker of high left ventricular end-diastolic pressure in patients with symptomatic left ventricular dysfunction. *Am. Heart J.* 1998; 135: 825-832.
17. Dao Q, Krishnaswamy P, Kazanegra R, Harrison A, Amirnovin R, Lenert L, Clopton P, Alberto J, Hlavin P, Maisel A. Utility of B-type natriuretic peptide in the diagnosis of congestive heart failure in an urgent-care setting. *J. Am. Coll. Cardiol.* 2001; 37: 379-385.
18. Mukoyama M, Nakao K, Saito Y, Ogawa Y, Hosoda K, Suga S, Shirakami G, Jougasaki M, Imura H. Increased human brain natriuretic peptide in congestive heart failure. *New Engl. J. Med.* 1990; 323: 757-758.

19. Sagnella GA. Measurement and significance of circulating natriuretic peptides in cardiovascular disease. *Clin. Science* 1998; 95: 519-529.
20. McDonagh TA, Robb SD, Murdoch DR, Morton JJ, Ford I, Morrison CE, Tunst 全部-Pedoe H, McMurray JJV, Dargie HJ. Biochemical detection of left-ventricular systolic dysfunction. *Lancet* 1998; 351: 9-13.
21. Mair J, Friedl W, Thomas S, Puschendorf B. Natriuretic Peptides in assessment of left-ventricular dysfunction. *Scand. J. Clin. Lab. Invest.* 1999; 59: 132-142.
22. Muders F, Kromer EP, Griese DP, Pfeifer M, Hense HW, Riegger GAJ, Elsner D. Evaluation of plasma natriuretic peptides as markers for left ventricular dysfunction. *Am. Heart J.* 1997; 134: 442-449.
23. Cowie MR, Struthers AD, Wood DA, Coats AJS, Thompson SG, Poole-Wilson PA, Sutton GC. Value of natriuretic peptides in assessment of patients with possible new heart failure in primary care. *Lancet* 1997; 350: 1347-1351.
24. Maisel AS, Krishnaswamy P, Nowak RM, McCord J, Hollander JE, Duc P, Omland T, Storrow AB, Abraham WT, Wu AH, Clopton P, Steg PG, Westheim A, Knudsen CW, Perez A, Kazanegra R, Herrmann HC, McCullough PA. Breathing Not Properly Multinational Study Investigators. Rapid measurement of B-type natriuretic peptide in the emergency diagnosis of heart failure. *N. Engl. J. Med.* 2002; 347: 161-167.
25. McCullough PA, Nowak RM, McCord J, Hollander JE, Herrmann HC, Steg PG, Duc P, Westheim A, Omland T, Knudsen CW, Storrow AB, Abraham WT, Lamba S, Wu AH, Perez A, Clopton P, Krishnaswamy P, Kazanegra R, Maisel AS. B-type natriuretic peptide and clinical judgment in emergency diagnosis of heart failure: analysis from Breathing Not Properly (BNP) Multinational Study. *Circulation* 2002; 106: 416-422.
26. Maisel AS, Koon J, Krishnaswamy P, Kazanegra R, Clopton P, Gardetto N, Morrissey R, Garcia A, Chiu A, De Maria A. Utility of B-natriuretic peptide as a rapid, point-of-care test for screening patients undergoing echocardiography to determine left ventricular dysfunction. *Am. Heart J.* 2001; 141: 367-374.
27. Lubien E, DeMaria A, Krishnaswamy P, Clopton P, Koon J, Kazanegra R, Gardetto N, Wanner E, Maisel AS. Utility of B-natriuretic peptide in detecting diastolic dysfunction. *Circulation* 2002; 105: 595-601.
28. Krishnaswamy P, Lubien E, Clopton P, Koon J, Kazanegra R, Wanner E, Gardetto N, Garcia A, DeMaria A, Maisel AS. Utility of B-natriuretic peptide in identifying patients with left ventricular systolic or diastolic dysfunction. *Am. J. Med.* 2001; 111: 274-279.
29. Omland T, Aakvaag A, Bonarjee VVS, Caidahl K, Lie RT, Nilsen DWT, Sundstord JA, Dickstein K. Plasma brain natriuretic peptide as an indicator of left ventricular systolic function and long-term survival after acute myocardial infarction. *Circulation* 1996; 93: 1963-1969.
30. Richards AM, Nicholls MG, Yandle TG, Ikram H, Espiner EA, Turner JG, Buttmore RC, Lainchbury JG, Elliott JM, Frampton C, Crozier IG, Smyth DW. Neuroendocrine prediction of left ventricular function and heart failure after acute myocardial infarction. *Heart* 1999; 81: 114-120.
31. Stein BC, Levin RI. Natriuretic peptides: physiology, therapeutic potential, and risk stratification in ischemic heart disease. *Am. Heart J.* 1998; 135: 914-923.
32. W 全部 en T, Landahl S, Hedner T, Nakao K, Saito Y. Brain natriuretic peptide predicts mortality in the elderly. *Heart* 1997; 77: 264-267.
33. Darbar D, Davidson NC, Gillespie N, Choy AMJ, Lang CC, Shyr Y, McNeill GP, Pringle TH, and Struthers AD. Diagnostic value of B-type natriuretic peptide concentrations in patients with acute myocardial infarction. *Am. J. Cardiol.* 1996; 78: 284-287.

34. Broome HE, Cembrowski GS, Kahn SN, Martin PL, Patrick CA. Implementation and use of a manual multi-rule quality control procedure. *Lab Med* 1985; 16: 533-537.
35. Westgard JO, Barry PL, Hunt MR, Groth T. A multi-rule Shewhart chart for quality control in clinical chemistry. *Clin Chem* 1981; 27: 493-501.
36. Koch DD, Ory 全部 JJ, Quam EF, Feldbruegger DH, et al. Selection of medic 全部 y useful QC procedures for individual tests done in a multitest analytical system. *Clin Chem* 1990; 36:230-233.
37. Mugan K, Carlson IH, Westgard JO. Planning QC procedures for immunoassays. *J Clin Immunoassay* 1994;17:216-222.
38. Approved Guidelines – Statistical Quality Control for Quantitative Measurements: Principles and Definitions, C24-A2. February 1999. National Committee for Clinical Standards.
39. Doust, J.A., Pietrzak, E., Dobson, A., and Glasziou, P., How well does B-type natriuretic peptide predict death and cardiac events in patients with heart failure: systematic review. *BMJ* 330:625-633, 2005.
40. Vrtovec, B., Delgado, R., Zewail, A., Thomas, C.D., Richartz, B.M., and Radovancevic, B., Prolonged QTc interval and high B-type natriuretic peptide levels together predict mortality in patients with advanced heart failure. *Circulation* 107:1764-1769, 2003.
41. Harrison, A., Morrison, L.K., Krishnaswamy, P., Kazanegra, R., Clopton, P., Dao, Q., Hlavin, P., and Maisel, A.S., B-type natriuretic peptide predicts future cardiac events in patients presenting to the emergency department with dyspnea. *Ann. Emerg. Med.* 39:131-138, 2002.
42. Cheng, V., Kazanegra, r., Garcia, A., Lenert, L., Krishnaswamy, P., Gardetto, N., Clopton, P., and Maisel, A., A rapid bedside test for B-type natriuretic peptide predicts treatment outcomes in patients admitted for decompensated heart failure: a pilot study. *J. Am. Coll. Cardiol.* 37:386-391, 2001.
43. Bettencourt, P., Ferreira, S., Azevedo, A., and Ferreira, A., Preliminary data on the potential usefulness of B-type natriuretic peptide levels in predicting outcome after hospital discharge in patients with heart failure. *Am. J. Med.* 2002 113:215-219, 2002.
44. Logeart, D., Thabut, G., Jourdain, P., Chavelas, C., Beyne, P., Beauvais, F., Bouvier, E., and Solal, A.C., Pre-discharge B-type natriuretic peptide assay for identifying patients at high risk of re-admission after decompensated heart failure. *J. Am. Coll. Cardiol.* 43:635-41, 2004.

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SDS

安全性数据表载于 [beckmancoulter.com/techdocs](http://beckmancoulter.com/techdocs)

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1. Cembrowski GS, Carey RN. Laboratory quality management: QC  $\hat{=}$  QA. ASCP Press, Chicago, IL, 1989.
2. Broome HE, Cembrowski GS, Kahn SN, Martin PL, Patrick CA. Implementation and use of a manual multi-rule quality control procedure. Lab Med 1985; 16: 533-537.
3. Westgard JO, Barry PL, Hunt MR, Groth T. A multi-rule Shewhart chart for quality control in clinical chemistry. Clin Chem 1981; 27: 493-501.
4. Koch DD, Oryall JJ, Quam EF, Feldbruegger DH, et al. Selection of medically useful QC procedures for individual tests done in a multitest analytical system. Clin Chem 1990; 36: 230-233.
5. Muga K, Carlson IH, Westgard JO. Planning QC procedures for immunoassays. J Clin Immunoassay 1994; 17: 216-222.
6. Approved Guidelines – Statistical Quality Control for Quantitative Measurements: Principles and Definitions, C24-A2. February 1999. National Committee for Clinical Standards.
7. Approved Guideline - Protection of Laboratory Workers From Occupationally Acquired Infections, M29-A4, 4th Edition, May 2014. Clinical and Laboratory Standards Institute.

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